

JHU NMR Experiment Quick Guide

1-Dimensional Experiments

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
Proton		proton		zg	R400/N400: 01m 50s R300: 02m 21s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)		FID Size (td, F2)
16	2	1	1		65536
Experiment Setup					
- Use getprosol to read in proper pulse parameters					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence		Experimental Time	
Proton-solv_supp		--		zgcppr		R400/N400: 10m 49s R300: 13m 39s	
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		1D Loop (td0)		FID Size (td, F2)	
128	2	1		1		65536	
Experiment Setup							
<div><div>-</div>Use getprosol to read in proper pulse parameters</div>							
<div><div>-</div>Move o1p to be on-resonance with the solvent peak</div>							

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
Deuterium		deuterium		zg2h	R400/N400: 04m 54s R300: 04m 02s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		1D Loop (td0)	FID Size (td, F2)
64	4	2		1	16384
Experiment Setup					
<div><div>-</div><div>This setup is through the lock channel, sensitivity will be lower compared to detection through the broad band channel.</div></div> <div><div>-</div><div>Open a <i>Proton</i> experiment and lock, tune, shim on the sample</div></div> <div><div>-</div><div>Once shimmed, read in the <i>Deuterium</i> parameters and perform rga (DO NOT reshim or tune)</div></div>					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
Boron11		boron11		zgpg30	R400/N400: 06m 06s R300: 07m 56s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)		FID Size (td, F2)
256	4	0.1	1		65536
Experiment Setup					
- Use getprosol to read in proper pulse parameters for decoupling					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
Carbon13-apt		c13apt		jmod	R400/N400: 14m 45s R300: 16m 42s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		1D Loop (td0)	FID Size (td, F2)
64	4	2		4	65536
Experiment Setup					
<div><div>-</div><div>Set cnst2 = J(XH), <i>J</i>-coupling between X and H.</div><div>-</div><div>If primarily aliphatic protons: cnst2=125.</div><div>-</div><div>If primarily aromatic protons: cnst2 = 150.</div><div>-</div><div>Default value: cnst2 = 140.</div><div>-</div><div>Set cnst11=1</div></div>					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
Carbon13-dec30		c13cpd		zgpg30	R400/N400: 14m 49s R300: 16m 46s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		1D Loop (td0)	FID Size (td, F2)
64	4	2		4	65536
Experiment Setup					
- Use getprosol to read in proper pulse parameters for decoupling					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence		Experimental Time	
Carbon13-dept90		-		Udeft_13C		R400/N400: 14m 57s R300: -----	
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		1D Loop (td0)		FID Size (td, F2)	
64	8	2		4		65536	
Experiment Setup							
<div>- Set cnst2 = J(XH), <i>J</i>-coupling between C and H.</div> <div>- Default value: cnst2 = 140.</div>							

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
Carbon13-dept135		-		Udeft_13C	R400/N400: 14m 57s R300: -----
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		1D Loop (td0)	FID Size (td, F2)
64	8	2		4	65536
Experiment Setup					
<div><div>-</div>Set cnst2 = J(XH), <i>J</i>-coupling between C and H.<div>-</div>Default value: cnst2 = 140.</div>					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
Carbon13-udeft		-		Udeft_13C	R400/N400 : 31m 02s R300 : 35m 22s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)	FID Size (td, F2)	
64	0	3.5	4	65536	
Experiment Setup					
<div><div>-</div><div>Do not change d1</div><div>-</div><div>d20 is an initial delay of the sequence which will only apply once. Should be set to be 5×<i>T</i>₁. Default value: d20 = 240s</div></div>					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time	
Nitrogen15		N15zg	zg	R400/N400: 06m 30s R300: 06m 39s	
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)		FID Size (td, F2)
32	4	10	1		32768
Experiment Setup					
- Use getprosol to read in proper pulse parameters					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time	
Nitrogen15-dec		N15cpd	zgig	R400/N400: 06m 30s R300: 06m 40s	
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)	FID Size (td, F2)	
32	4	10	1	32768	
Experiment Setup					
- Use getprosol to read in proper pulse parameters for decoupling					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
Nitrogen15-inept		N15inept		ineptrd	R400/N400: 08m 40s R300: 08m 53s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)	FID Size (td, F2)	
32	16	10	1	32768	
Experiment Setup					
<div><div>-</div><div>Set cnst2 = J(XH), <i>J</i>-coupling between N and H.</div><div>-</div><div>Default value: cnst2 = 90.</div><div>-</div><div>Set cnst11=6</div></div>					

Current <i>rpar</i> Name	Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time	
Silicon29-dept	si29dept	dept45	R400/N400 : 51m 02s R300 : 53m 01s	
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)	FID Size (td, F2)
128	4	5	4	65536
Experiment Setup				
<div><div>-</div>Set cnst2 = J(XH), <i>J</i>-coupling between Si and H.<div>-</div>Default value: cnst2 = 7</div>				

Current <i>rpar</i> Name		Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time
Silicon29-udefit		-	Udefit_29Si	R400/N400: 49m 52s R300: -----
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)	FID Size (td, F2)
128	0	3.5	4	65536
Experiment Setup				
<div><div>-</div><div>Do not change d1</div></div> <div><div>-</div><div>d20 is an initial delay of the sequence which will only apply once. Should be set to be 5×T₁. Default value: d20 = 90s</div></div>				

Current <i>rpar</i> Name		Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time	
Phosphorus31-dec30		p31cpd	zgpg30	R400/N400: 05m 29s R300: 05m 31s	
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)		FID Size (td, F2)
128	4	2	1		65536
Experiment Setup					
- Use getprosol to read in proper pulse parameters for decoupling					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time	
Selenium77		s77_1D	zg	R400/N400: 05m 56s R300: -----	
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)		FID Size (td, F2)
1024	4	0.1	1		16384
Experiment Setup					
- Use getprosol to read in proper pulse parameters					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time	
Selenium77-dec		s77cpd	zgpg	R400/N400: 06m 27s R300: -----	
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)		FID Size (td, F2)
1024	4	0.1	1		16384
Experiment Setup					
- Use getprosol to read in proper pulse parameters for decoupling					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time	
Fluorine19		F19zg	zgflqn	R400/N400: ----- R300: 04m 11s	
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)		FID Size (td, F2)
64	4	1	1		131072
Experiment Setup					
- Use getprosol to read in proper pulse parameters					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time	
Fluorine19-dec		F19cpd	zgfhigqn	R400/N400: ----- R300: 04m 11s	
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)		FID Size (td, F2)
64	4	1	1		131072
Experiment Setup					
- Use getprosol to read in proper pulse parameters and decoupling parameters					

1-Dimensional Correlation Experiments

Current <i>rpar</i> Name		Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time
NOESY-sel_1D		selNOE_1D/1D_gradNOE	selnoggp	R400/N400: 10m 58s R300: 13m 46s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)	FID Size (td, F2)
32	4	2	4	32768
Experiment Setup				
<ul style="list-style-type: none"> - d8 = mixing time. - Small molecules: $d8 \approx T_1$ - Large molecules: $d8 \approx 0.5T_1$ - Default d8 = 800ms - spoffs2 = offset frequency for selective shaped pulse - set to the frequency difference between the center (o1) and the location of the desired peak to be analyzed. 				

Current <i>rpar</i> Name		Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time
NOESY-pr_1D		solvent_suppression	noesypr1d	R400/N400: 10m 45s R300: 13m 45s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)	FID Size (td, F2)
32	4	2	4	32768
Experiment Setup				
<ul style="list-style-type: none"> - d8 = mixing time. - Small molecules: $d8 \approx T_1$ - Large molecules: $d8 \approx 0.5T_1$ - Default d8 = 300ms 				

Current <i>rpar</i> Name		Previous <i>rpar</i> Name	Pulse Sequence	Experimental Time
TOCSY-sel_1D		selTOCSY_1D	selmlgp	R400/N400: 02m 51s R300: 03m 18s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	1D Loop (td0)	FID Size (td, F2)
32	4	2	1	32768
Experiment Setup				
<ul style="list-style-type: none"> - d9 = mixing time. - Default: d9 = 60 ms - spoffs2 = offset frequency for selective shaped pulse - set to the frequency difference between the center (o1) and the location of the desired peak to be analyzed. 				

2-Dimensional Experiments

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
COSY-dqf		-		cosydfph	R400/N400: 07m 36s R300: 08m 05s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	FID Size (td, F2)		2D Slices (td, F1)
1	16	1.5	2048		256
Experiment Setup					
- Use getprosol to read in proper pulse parameters					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
COSY-grad		cosygpsw		cosygpqf	R400/N400: 07m 39s R300: 08m 10s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		FID Size (td, F2)	2D Slices (td, F1)
1	16	1.5		2048	256
Experiment Setup					
- Use getprosol to read in proper pulse parameters					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
COSY-HF		HFcorrelation		hfcoqfq	R400/N400: ----- R300: 33m 25s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	FID Size (td, F2)		2D Slices (td, F1)
16	16	1.0	2048		256
Experiment Setup					
- Use getprosol to read in proper pulse parameters					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
NOESY		NOESY		noesyph	R400/N400: 1h 25m 05s R300: 1h28m 50s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		FID Size (td, F2)	2D Slices (td, F1)
8	16	2		2048	256
Experiment Setup					
<div>- d8 = mixing time.</div> <div>- Small molecules: $d8 \approx T_1$</div> <div>- Large molecules: $d8 \approx 0.5T_1$</div> <div>- Default d8 = 300ms</div>					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
NOESY-HF		HFnoesy		hoesyfqph.cdm	R400/N400: ----- R300: 33m 25s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		FID Size (td, F2)	2D Slices (td, F1)
16	16	1.0		4096	256
Experiment Setup					
<div>- d8 = mixing time.</div> <div>- Default d8 = 400ms</div>					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
NOESY-pr		noesypr		noesyphpr	R400/N400 : 1h 26m 05s R300 : 1h 29m 50s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	FID Size (td, F2)		2D Slices (td, F1)
8	16	2	2048		256
Experiment Setup					
<div><div>-</div><div>d8 = mixing time.</div><div>-</div><div>Small molecules: $d8 \approx T_1$</div><div>-</div><div>Large molecules: $d8 \approx 0.5T_1$</div><div>-</div><div>Default d8 = 300ms</div></div>					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
NOESY-grad		-		noesygpqh	R400/N400 : 1h 06m 54s R300 : 1h 28m 16s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	FID Size (td, F2)	2D Slices (td, F1)	
16	16	0.5	2048	256	
Experiment Setup					
<div><div>-</div><div>d8 = mixing time.</div><div>-</div><div>Small molecules: $d8 \approx T_1$</div><div>-</div><div>Large molecules: $d8 \approx 0.5T_1$</div><div>-</div><div>Default d8 = 300ms</div></div>					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence		Experimental Time	
TOCSY		TOCSY		mlevph		R400/N400: 1h 16m 42s R300: 1h 20m 36s	
Scans (ns)	Dummy Scans (ds)		Recycle Delay (d1)		FID Size (td, F2)		2D Slices (td, F1)
8	16		2		2048		256
Experiment Setup							
<ul style="list-style-type: none">- d9 = mixing time, longer the delay further the coupling distance<ul style="list-style-type: none">- Default: d9 = 60 ms- spoffs2 = offset frequency for selective shaped pulse<ul style="list-style-type: none">- Set to the frequency difference between the center (o1) and the location of the desired peak to be analyzed.							

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
HSQC-ed		edited-hsqc		hsqcedetgp	R400/N400: 15m 05s R300: 15m 42s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)	FID Size (td, F2)	2D Slices (td, F1)	
2	16	1.5	2048	256	
Experiment Setup					
<div><div>-</div><div>cnst2 = J(XH)</div><div>-</div><div>Default: cnst2 = 145 Hz</div><div>-</div><div>d21 = multiplicity selection by setting to 1/(2J(XH))</div></div>					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
HMBC-grad		gradientHMBC		hmbcgplpndqf	R400/N400: 1h 07m 18s R300: 1h 11m 51s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		FID Size (td, F2)	2D Slices (td, F1)
16	16	1.5		4096	128
Experiment Setup					
<div><div>-</div>Set cnst2 and cnst13 to the short and long range X-H J-couplings</div> <div><div>-</div>cnst2 = J(XH)<div><div>-</div>Default: cnst2 = 145 Hz</div></div> <div><div>-</div>cnst13 = J(XH) long range couplings<div><div>-</div>Default: cnst13 = 10 Hz</div></div>					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
HMBC-j_ed		j-HMBC		hmbcetgpl3nd	R400/N400: 1h 25m 11s R300: 1h 29m 17s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		FID Size (td, F2)	2D Slices (td, F1)
16	16	2		4096	128
Experiment Setup					
<div><div>-</div><div>Set cnst6 and cnst7 to the minimum and maximum X-H 1-bond J-couplings</div><div>-</div><div>cnst6/cnst7 = 1J(XH)</div><div>-</div><div>Default: cnst6 = 110 Hz</div><div>Cnst7 = 170 Hz</div><div>-</div><div>cnst13 = J(XH) long range couplings</div><div>-</div><div>Default: cnst13 = 10 Hz</div></div>					

Current <i>rpar</i> Name		Previous <i>rpar</i> Name		Pulse Sequence	Experimental Time
H1-Kinetics		-		zg2dcppr	R400/N400 : 10h 40m 03s R300 : 10h 40m 04s
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		FID Size (td, F2)	2D Slices (td, F1)
8	0	10		65536	128
Experiment Setup					
<div>- d20 = delay between experiments</div> <div>- Default : d20 = 300s</div>					

Previous <i>rpar</i> Name		Current <i>rpar</i> Name		Pulse Sequence	Experimental Time
DOSY-ledbpgp2s		-		ledbpgp2s	R400: varies R300: -----
Scans (ns)	Dummy Scans (ds)	Recycle Delay (d1)		FID Size (td, F2)	2D Slices (td, F1)
8	0	**		16384	32
Experiment Setup					
<ul style="list-style-type: none">- Setup and acquire a standard proton experiment.<ul style="list-style-type: none">- Set the sw and center your data in the spectrum so that there is approximately 10% of clear baseline (no signal) on both sides signal region.- Measure the T_1 of your sample (see T_1 measurement notes) and set the d1 = 5×T_1 with respect to the longest relaxation measurement.- Read in the DOSY-ledbpgp2s experiment- Spin your sample at 20 Hz- In the command line enter the following to start the experiment: xau dosy 5 80 32 l y y					